
AQUA PULSE

- ES.** Consejos y advertencias de seguridad
- EN.** Safety advice and warnings
- FR.** Conseils et avertissements de sécurité
- PT.** Conselhos e advertências de segurança
- IT.** Consigli e avvertimenti di sicurezza
- CA.** Consells i advertències de seguretat
- DE.** Ratschläge und sicherheitshinweise
- RO.** Recomandări siguranță și avertizări
- PL.** Wskazówki i ostrzeżenia bezpieczeństwa
- BG.** Съвети и предупреждения за безопасност
- AR.** نصائح وتحذيرات تتعلق بالسلامة

solac

ES

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DENTAL IRRIGATOR AQUA PULSE

Dear customer,

Many thanks for choosing to purchase a SOLAC brand product.

Thanks to its technology, design and operation and the fact that it exceeds the strictest quality standards, a fully satisfactory use and long product life can be assured.

SAFETY ADVICE AND WARNINGS

Read these
instructions

carefully before

switching on the appliance and keep them for future reference. Failure to follow and observe these instructions could lead to an accident.

- Oral irrigators can be used by children aged from 8 years and above and persons with reduced physical, sensory

or mental capabilities or lack of experience and knowledge under supervision.

- This oral irrigator shall be used only with cold or warm water (below 40°C).
- Children shall not play with the appliance.
- This appliance is for household use only, not professional or industrial use.

This appliance has a detachable supply unit.
(ZD006M050100EU).

- This appliance is only to be used with the power supply unit provided with the appliance.
- This appliance contains batteries that are non-replaceable. When the

battery is at end of life, the appliance shall be properly disposed of.

- **WARNING:** For the purposes of recharging the battery, only use the detachable supply unit provided with this appliance.
- Ensure that the voltage indicated on the rating label matches the mains voltage before plugging in the appliance.
- The appliance's plug must fit into the mains socket properly. Do not alter the plug. Do not use plug adaptors.
- Do not use the appliance if the detachable supply unit is damaged.
- If the appliance enclosure breaks, immediately disconnect the appliance from the mains to prevent the possibility of an electric shock.
- Do not touch the plug with wet hands.

SERVICE

- Any misuse or failure to follow the instructions for use renders the guarantee and the manufacturer's liability null and void.

ANOMALIES AND REPAIR

- Take the appliance to an authorised technical support service if problems arise. Do not try to dismantle or repair

without assistance, as this may be dangerous.

FOR EU PRODUCT VERSIONS AND/OR IN THE CASE THAT IT IS REQUESTED IN YOUR COUNTRY: ECOLOGY AND RECYCLABILITY OF THE PRODUCT

- The materials of which the packaging of this appliance consists are included in a collection, classification and recycling system. Should you wish to dispose of them, use the appropriate public containers for each type of material.
- The product does not contain concentrations of substances that could be considered harmful to the environment.

This symbol means that in case you wish to dispose of the product once its working life has ended, take it to an authorised waste agent for the selective collection of waste electrical and electronic equipment (WEEE).

- This appliance complies with Directive 2014/35/EU on Low Voltage, Directive 2014/30/EU on Electromagnetic Compatibility, Directive 2011/65/EU on the restrictions of the use of certain hazardous substances in electrical and electronic equipment and Directive 2009/125/EC on the ecodesign requirements for energy-related products.
- The following information is regarding regulation 2019/1782 on ecodesign

requirements for external power
supplies:

Trademark and address	Solac Small Electro, S.L. Avda. Barcelona, s/n, 25790 Oliana, Lleida (Spain)
Model	ZD006M050100EU
Input voltage	100-240 V
Input AC frequency	50/60 Hz
Output voltage	5,0 V
Output current	1,0 A
Output power	5,0 W
Average active efficiency	74,00%
Efficiency at low load (10 %)	N.A.
No-load power consumption	0,09 W

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ENGLISH

WARRANTY AND TECHNICAL ASSISTANCE

- This product enjoys the recognition and protection of the legal guarantee in accordance with current legislation. To enforce your rights or interests you must go to any of our official technical assistance services.
- You can find the closest one by accessing the following web link: <http://solac.com/>
- You can also request related information by contacting us.

- You can download this instruction manual and its updates at <http://solac.com/>

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